

Wheeled Coach Stock Series



CitiMedic Plus

Key Features 2100 Lb. Payload

Printed Circuit Board
LED self-diagnostics.
Eliminates more than 300 wire crimps and connections to minimize failure areas. Lifetime electrical warranty standard.

More Room
Each CitiMedic comes with a spacious interior and ergonomically designed working environment for emergency care.

Corner Cap Lights
Wheeled Coach's corner caps provide structural closure of roof rail tangent points, welded at each corner of the modular body. Also available with optional strobe or LED package.

AVW Grille Assembly
Patented design maximizes air flow for cooling systems and protects your speakers and warning lights.

CHASSIS

- Ford E-350 Dual Rear Wheel Cutaway, 138" WB, 11,500# GVWR
- 6.0 L Turbo-Diesel (445 CID), Ambulance Prep Package
- 4600# Front Axle, 7800# Rear Axle
- Electronic 4-speed Automatic Overdrive Wide Ratio Transmission
- Dual OEM Alternators
- Throttle Idler
- Dual OEM Batteries, 1700 CCA Total
- 4.10 Limited-Slip Axle Ratio
- Auxiliary Transmission Oil Cooler
- 684 Handling Package
- 57L Auxiliary Heater and A/C Package
- Tilt and Cruise
- Intermittent Wipers
- AM/FM Electronic Stereo with Clock
- High Back Velour Seats with Armrests (XL Pkg)
- Power Windows and Door Locks (XL Pkg)
- Carpet on Cab Floor (XL (Pkg))
- 37 Gallon Fuel Tank
- Driver Air Bag

PATIENT MODULE

Module:

- 144" X 90"

Headroom:

- 68" Interior with Crowned Roof to Inhibit Water Pooling

Exterior Compartments:

- Right Rear Compartment, Spare Tire Storage
- Left Front Compartment, Oxygen Storage with Status Board
- Vertical Backboard Storage, Left Rear Compartment

Windows:

- Curbside Door, Sliding, Tinted Glass with Screen

- Rear Doors, Fixed, Tinted

Rear Bumper:

- Non-Skid Expanded Metal with Center Flip-up and Under Body Skids with Integral Tow Eyes

Fender Flare:

- Black Rubber

Skirtrails:

- Lower, Diamond Plate

Kickplate:

- Diamond Plate with Recessed Tag Area

Paint Design: (Optional)

Insulation:

- Fiberglass in Module and Ceiling, Solid Block Foam in Doors

EMERGENCY LIGHTS

Light Bar: (Optional)

Flasher:

- Dual Mode Solid State

Grille Lights:

- (2) 7 x 3 Halogen, Red Lens

Intersection Lights:

- (2) 7 x 3 Halogen, Red Lens

Warning Lights:

- (8) 9 x 7 Halogen, Red Lens (2) Front, (2) Rear, (2) Each Side
- (1) 7 x 3 Halogen, Amber, Center Rear
- (1) 9 x 7 Halogen, Clear, Center Front

Load Lights:

- (2) Optiscene, 9 x 7, 26° Angle on Rear

Scene Lights:

- (4) Optiscene, 9 x 7, 13° Angle (2) Each Side

Tail Lights:

- (6) Round Sealed Lights Installed in Rear Diamond Plate Kick Panel

Siren:

- Whelen WS-295HFS, Located in Front Switch Console

Siren Speakers:

- Wheeled Coach Grille AVW System with (2) 100 Watt Drivers and the (2) Grille Lights

Back-Up Alarm:

- Installed Under Rear of Vehicle. Cut-Off Switch (Optional)

DOT Clearance Lights:

- Exclusive Integrated Corner-Cap Lights and (3) Amber Front and (3) Red Rear

Spotlight:

- (1) Handheld, 100,000 CP

Patient Domes:

- (5) Weldon Model 8045 Halogen Lights Mounted in Headliner, (3) Over Cot, (3) Over S/B

OXYGEN/VACUUM

Oxygen Outlets:

- (2) Ohio Style Fitting in A/A

Vacuum Pump:

- Thomas E-12V, Mounted Under Attendant Seat

CABINETS & HARDWARE

Radius Corner Wood Cabinets & Seamless Upholstery

Right Front ALS Cabinet:

- Designed for Storage of Monitor and Multiple Medical Cases

Squadbench:

- 72" Bench with (1) Seat Cushion, and (3) Individual Backrests, Sleeved Seat Belts for (3) Sitting Positions or (1) Secondary Patient

Infectious Waste Cabinet:

- Slide Out Trash Receptacle
- Sharps Container Storage

Basewall Cabinets:

- (2) Full Depth Cabinets Across the Top with Sliding Plexiglass Windows and (4) Full Size Cabinets with (1) Adjustable Shelf in Each
- Stair Chair Storage at Rear of Wall, Near Rear Doors

Attendant Seat Base:

- With Venting for Mounting of Electrical Components

Patient Module Liner:

- One Piece Fiberglass

Grabrail:

- 1" Diameter Handicap Style with Rounded Ends

Assist Handle:

- (1) Each Entry Door

IV Holders:

- (2), (1) Over Cot and (1) Over Squadbench Chest Areas

Attendant Seat:

- Cushion on Wood Cabinet

Plexiglass Doors:

- Sliding Doors with Full-Length Aluminum Handles

ELECTRICAL

Heat-Air Conditioning System:

- Filtered Air System
- Fully Automatic with Thermostat

Circuit Breakers:

- Auto Resetting, Pole Style on a Printed Circuit Board

Shoreline:

- STD 15 AMP GFI Protected Circuit, Inlet Mounted by Driver's Door on Front of Module

- (1) Outlet Mounted in A/A

- (1) Outlet Mounted on Curbside Wall, at Head of Squadbench

12 Volt Outlets:

- Power Point, (2) in A/A

STAR OF LIFE DECALS

(Optional)

